
MANUAL

Incubator P S1 compact

This operating manual is an integral part of this product. It contains important information on installing and using the product. Please keep this in mind, especially when passing this product on to third parties. For this reason, keep the operating instructions for future reference!

Design and specifications are subject to change without notice. The manual is not covered by an update service.

Unless expressly authorized, forwarding and duplication of this document, and the utilization and communication of its contents are not permitted. Violations will entail an obligation to pay compensation.

All rights reserved in the event of granting of patents or registration of a utility model.

All product names mentioned herein may be the trademarks or registered trademarks of their respective companies and "TM" and "®" are not mentioned in each case in this manual.

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© 2018

Contents

Description	1
Safety instructions and hazard warnings	2
Usage in accordance with intended purpose	3
Scope of delivery	4
Operating controls	4
Assembly	5
Operation	7
Maintenance	8
Troubleshooting	8
Specifications	9

Description

- The INCUBATOR P S1 COMPACT is designed for the stabilization of *in vitro* conditions for cell- and tissue cultures during microscopic examination.
- The incubator enables a homogenous heat, CO₂ and O₂ distribution in combination with a Heating Insert.
- The heated glass insert warms up the incubation chamber from above. This avoids the condensation of water on the cover of the cell cultivation vessel.
- The heated glass insert of the incubator is permeable to 90% of light in the visible wavelength range.
- The INCUBATOR P S1 COMPACT is also suitable for high-resolution microscopy. It is designed for the LD-condensors 0.35 and 0.55.
- For electrical power supply and the control of the temperature inside the incubator, the TEMPMODULE S has to be used.
- For the control of the CO₂ concentration, the CO₂ MODULE S has to be used.
- For the control of the O₂ concentration, the O₂ MODULE S has to be used.

Safety instructions and hazard warnings

= **Notes for user safety and to avoid damage to the device. Read the instructions for use.**

= *Advice for a trouble-free operation of the device.*

= *Further information for systems optimization.*

Due to the following conditions the glass can be hot:

- during heat-up phase
- ambient temperature below $\sim 19^{\circ}\text{C}$
- during cold air motion above the incubator

- Do not feed the incubator with combustible gases, vapors or liquids.
- Do not autoclave the incubator.
- If the incubator is used with local gassing solutions, a sufficient room ventilation has to be ensured to prevent the accumulation of a harmful CO_2 -concentration (starting from approx. 0.5 Vol-%) or the accumulation of a harmful N_2 -concentration (reducing the oxygen concentration).
- If you have reasons to assume that a safe operation is no longer possible, immediately take the device out of operation and secure it against unintentional use. Reasons to assume that safe operation is no longer possible include:
 - the device shows signs of visible damage,
 - the device does not work anymore,
 - the device has been stored for long periods under unfavourable conditions, or
 - has been subjected to considerable stress in transit.
- You should also follow the additional safety instructions in each chapter of this manual as well as in the manuals of the connected devices.

Usage in accordance with intended purpose

- The product is designed for laboratory use by qualified personnel only.
It is no medical product according to the IVD Directive 98/79/EG or the MD Directive 93/42/EEG. Therefore, this product is not allowed to be used for medical analysis, diagnostic or treatment purposes on humans within the scope of the IVD or MD Directive.
- The INCUBATOR P S1 COMPACT offers a homogenously heated incubation room with a defined atmosphere in combination with following inserts:

PeCon Article-No.	Insert
130-800 031	Heating Insert P S1 compact
131-800 032	Heating Insert P Lab-Tek™ S1 compact

- For power supply and the control of temperature, the CO₂ concentration and the O₂ concentration the INCUBATOR P S1 COMPACT has to be used with following devices:

PeCon Article-No.	Device	Parameter
800-800 000	TempModule S	°C (temperature)
810-800 001	CO ₂ Module S	CO ₂ (carbon dioxide)
820-800 002	O ₂ Module S	O ₂ (oxygen)
850-800 011	Heating Device Humidity S1	Humidity

- The heating power is designed for a stable heating up to 30°C above ambient temperature.
- The incubator is designed for increased CO₂ and decreased O₂ values.
- This product is designed for indoor use. Operation is not permitted under unfavourable ambient conditions. The following are unfavourable ambient conditions:
 - dampness or humidity
 - combustible
 - gases, vapours or solvents


- Use other than the described above can damage the product and other components used with it. It may additionally involve other risks. Do not change or modify any part of the product. The safety instructions must be observed without fail.**
- Any deviations from the application described here will forfeit all claims against warranty and the manufacturer's liability in general and particularly in the event of a defect. Furthermore, the manufacturer is not liable for consequential damage unless intention or gross negligence can be proved.**


Scope of delivery

		Qty.	Article-No.
	Incubator P S1 compact	1	#160-800 094
	Operating manual	1	--

Operating controls

- (1) Gas mixture inlet
- (2) Openings for cables, perfusion tubes and channel closures of different inserts
- (3) Aluminium base frame
- (4) Heatable glass
- (5) Electric connecting cable


- (6)** 8-pin connector
- (7)** Openings for gas mixture output
- (8)** Temperature sensor Pt100

For instructions within this manual, we refer again to the numbers of this legend **(x)**.

Assembly

Connection diagrams

Application for temperature control, CO₂ and O₂ control with the CO₂ and O₂ MODULE S:

Insertion of the insert into the stage

1. Put the Heating Insert (see section "Usage in accordance with intended purpose") into the cut-out (160 x 110 mm) of the scanning stage or mechanical stage.
Make sure that the spring in the microscope stage clamps the insert.
2. For a constant focus, adjust the insert with the four adjusting screws in the edges of the insert by using the supplied Allen key.

Attaching of the incubator

Put the INCUBATOR P S1 COMPACT onto the Heating Insert from above.

Electrical connection

1. Connect the electric cable of the Heating Insert:
 - Plug the 8-pin connector into the channel 1 of the TEMPMODULE S.
2. Connect the electric connecting cable **(5)** of the INCUBATOR P S1 COMPACT:
 - Plug the 8-pin connector **(6)** into the channel 2 of the TEMPMODULE S.

CO₂ control with the CO₂-Module S

The INCUBATOR P S1 COMPACT has to be connected to the CO₂ MODULE S and creates the precondition for an optimal pH-value in the culture medium.

Since CO₂ is heavier than the ambient air, it primarily escapes through leaks in the bottom area of the incubation system. Therefore, it is important that the cell cultivation vessels being used cover and seal up the observation opening of the baseplate completely.

1. Connect the CO₂ MODULE S to the CO₂ supply and the mains supply (see manual CO₂ MODULE S).
2. Connect the INCUBATOR P S1 COMPACT to the Humidification bottle (filled with aqua dest.) that is supplied with the CO₂ MODULE S:
 - Plug the transparent Tygon® tube (supplied with the CO₂ MODULE S) on the hose nipple (marked with "incubator") at the Humidification bottle.
 - Attach the other end of the transparent Tygon® tube to the gas mixture inlet **(1)** of the INCUBATOR P S1 COMPACT.

If you use the Heating Device Humidity S1, use the insulated tube instead of the transparent tube. The insulation largely avoids the condensation of water on the tube wall. The insulated tube is provided with the Heating Device Humidity S1.

3. Connect the Humidification bottle to the CO₂ MODULE S:
 - Plug the blue plastic tube (supplied with the CO₂ MODULE S) into the blue plug connector of the Humidification bottle.
 - Attach the other end of the blue plastic tube to the plug connector on the rear side of the CO₂ MODULE S (see manual CO₂ MODULE S).

Reduction of O₂ with the O₂ Module S

Arrange the modules according to the guideline "Small gas feeding". You find this guideline in the manual that is provided with the modules.

Small gas feeding (internal circulation):

Link the connector of the CO₂ MODULE S to the connector at the supplied 500 ml Humidification bottle (replaces the 250 ml Humidification bottle of the CO₂ MODULE S) with the 2m blue PU-tube (PUN-4). Link the Humidification bottle to the small incubator with the 2m clear HAT-tube.

To avoid the occurrence of a temperature gradient in the cell cultivation medium, mount the OBJECTIVE HEATER S1 onto the objective.

Operation

Switch-on

Every heated component has specific properties regarding their heating and cooling behavior. This has to be considered for a stable control of the temperature of this component. To obtain the best adaptation of the loop control to the heated component, the control parameters for channel 1 and channel 2 of the temperature control device have to be configured correctly.

Application with the TempModule S:

- Switch on TEMPMODULE S.
- With the software MTB2004 (part of software package Zeiss AxioVision 4.6 or higher) configure channel 1 with "**Device: H Insert P**", if a Heating Insert is used.
- Configure the channel 2 to "**Device: Inc PM**".
- With the software AxioVision 4.6 set the nominal temperature and activate the heating. Depending on the microscope used, the configuration can be done with the PC or TFT-panel. For detailed information see the manual of the software or TEMPMODULE S.
- Choose the appropriate parameters on the user interface. For proposals, see the manual for the TEMPMODULE S.
- During heat-up close the observation area with the red insert (see manual for inserts).
- If the nominal values are reached (duration approx. 5 min), the incubation system is ready for operation and the cell cultivation chamber can be inserted.
- If water condensates at the cover of the cell cultivation vessel, it can be removed by increasing the air temperature (increase the nominal value of the incubator).

For detailed information see the manual of TEMPMODULE S.

Maintenance

- The incubator is service free, apart from occasionally cleaning.
- Clean the incubator with a soft, lint-free, antistatic and light moistened cloth.
- For the cleaning of the heatable glass parts, we recommend water with mild detergents or alcohol (ethanol, 70%).
- For the cleaning of the plastic body, we recommend acrylic glass cleaners or water with mild detergents.

Do not use any abrasive or chemical detergents or detergents containing solvents.
Do not use cleaners with alcohol for cleaning the plastic body. This can lead to stress cracks.
The incubator cannot be autoclaved.

Troubleshooting

By purchasing this device, you have acquired a product that is reliable and operationally safe. However, problems and malfunctions may occur.

Periodically check the technical safety of the device, e.g. check for damage to the housing etc.

Do not fail to follow the safety instructions.

The following table shows possible problems, their causes and solutions.

Problem:	Cause / Solution:
1. <i>The INCUBATOR P S1 COMPACT does not heat.</i>	a) The TEMPMODULE S is not connected to mains power or is not switched on. Check the setup of the TEMPMODULE S (see manual TEMPMODULE S). b) The main fuse of the TEMPMODULE S is blown. Check the main fuse of the device and replace it if necessary (see manual TEMPMODULE S). c) The incubator is not connected to the TEMPMODULE S. Connect the INCUBATOR P S1 COMPACT with the TEMPMODULE S (see section "Assembly") d) The channels of TEMPMODULE S are not activated. Activate both channels of the TEMPMODULE S (see manual TEMPMODULE S).

Problem:	Cause / Solution:
2. <i>The INCUBATOR P S1 COMPACT does not reach the nominal temperature.</i>	a) It takes approx. 10 min to reach a working temperature of 37°C. b) The channel 1 is not configured with "H Insert P" (see "Operation"). Check the configuration on the TEMPMODULE S (see manual TEMPMODULE S). c) The channel 2 of the TEMPMODULE S is not configured with "Inc PM" (see "Operation"). Check the configuration on the TEMPMODULE S (see manual TEMPMODULE S).
3. <i>The cells already die at an early point of time or show unusual behaviour.</i>	The gas mixture output of the CO₂ MODULE S or the O₂ MODULE S is blocked. Clean it with a higher air pressure (max. 2.5 bar/36 psi).

Any other repair must always be carried out by qualified experts familiar with the hazards involved and with the relevant regulations. Unauthorized modifications and repairs of or inside the device will render the warranty null and void.

In case of further questions consult the technical information service of Carl Zeiss Micro Imaging GmbH:

e-mail: mikro.service@Zeiss.de
web: <http://www.zeiss.de/mikro>

Beyond that our technical information service would be glad to answer your questions:

e-mail: support@pecon.biz
web: <http://support.pecon.biz>

Specifications

Material	Acrylic glass (PMMA) black, glass
DC operating voltage	24 V dc protective low voltage
Power consumption	max. 0.5 A
Temperature sensor	Pt100
Control range	Ambient temperature up to 30°C above ambient
Dimensions (WxHxD)	116 x 12.5 x 190 (mm)
Observation area (LxW)	120 x 77 mm
Weight	approx. 0.3 kg

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We hereby declare to the subject matter described below, the agreement with the

- Directive 2006/95/EC of the European Parliament and of the Council of 12 December 2006 on the harmonisation of the laws of Member States relating to electrical equipment designed for use within certain voltage limits,
- Directive 2004/108/EC of the European Parliament and of the Council of 15 December 2004 on the approximation of the laws of the Member States relating to electromagnetic compatibility,
- and the Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances (RoHS) in electrical and electronic equipment.

If modification to the product not authorized by us, this declaration loses its validity.

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The declaration certifies the compliance with the directives and the law.
Conditions of guarantee and liability are dealt within our General Conditions of Sale.

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PeCon GmbH