

Heating unit:

Avoid excessive Induction noise, Static electricity, magnetic fields

• Parts: base, lid, objective

• Senses ambient temperature and can monitor culture temperature with thermocouple ^{→ around sample}

• Connectivity with PC

↳ Connectivity

with OKA-Touch (touch panel to control parameters)

↳ Can connect to separate units such as Bold Line Gas Controller

• Chamber feedback mode monitors base and lid temp and combines with pre-defined/calibrated ambient temp to get sample to closest temp. (e.g. knows at ambient $\approx 23^{\circ}\text{C}$; Need base at 38°)

• Allows for temperature cycling

Technical Specs: • Chamber Temp range: $25 - 60^{\circ}\text{C}$

Sample Temp range: Ambient $+3 - 60^{\circ}\text{C}$

Resolution: 0.1°C

Accuracy: 0.1°C in sample

0.2°C in chamber midp

Rise time 30 mins

• operate $0^{\circ}\text{C} - 55^{\circ}\text{C}$

• Store $-5^{\circ}\text{C} - 60^{\circ}\text{C}$

• operate humidity $0 - 20^{\circ}$

• weight: 1.6 Kg

CO₂/O₂ unit

multiple variants for different CO₂/O₂ control ranges

can use pressurized CO₂/O₂/N₂ for full control of CO₂ and O₂

alternatively can use CO₂ + air pump for just CO₂ control

- Requires PC or optional control panel for operation
- supports gas cycling
- Interfaces with T unit to allow communication between devices

